

Place the operating software (OtoTestUSB.exe) in a directory that you choose or create one yourself (e.g. C:\usbleakagetester) and make a shortcut.  
If, after starting the app, a security alert comes up you can avoid this in the future by unchecking the box.

The files logo.jpg (.png) and certifcate.pdf should be placed in the same directory as the operating software.

For this app it is no longer possible to connect two or more leakage testers to the same computer.

The maximum size of the logo should not exceed 320x200 pixels. It is best to choose a white background and the filename.ext must be "logo.jpg" or "logo.png"; if these logo-files are not present at all the Cursor logo (EARS) will be displayed.

Your own certificate design must be marked as "certificate.pdf", otherwise the alert "File not found" will appear.